

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human FVIIa Mutant Version (Vatreptacog alfa) Antibody (4F5)
PRODUCT REFERENCE	ARO-A17401
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	FVIIa Mutant Version (FVIIaDVQ)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-FVIIa Mutant Version (FVIIaDVQ), Anti-FVIIaDVQ, Anti-FVIIa Mutant, Vatreptacog alfa