



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |  |
|--------------------------|--|
| <b>PRODUCT NAME</b>      | Anti-Human NMU/Neuromedin-U Antibody (ABS0155) |
| <b>PRODUCT REFERENCE</b> | ARO-A17615                                     |
| <b>CLONALITY</b>         | Monoclonal Antibody                            |
| <b>ISOTYPE</b>           | IgG2a, kappa                                   |
| <b>TARGET</b>            | NMU  |

### PRODUCT INFORMATION

---

|                                |   |
|--------------------------------|---|
| <b>ORIGIN SPECIES</b>          | <i>Mus musculus (Mouse)</i>   |
| <b>FORM</b>                    | Liquid  |
| <b>BUFFER</b>                  |   |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term.<br>Addition of 20 to 40 % glycerol may improve cryoprotection.<br>Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.<br>Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | Anti-NMU, Anti-Neuromedin-U   |

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**