

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human TOR1A/Torsin-1A/DQ2 Antibody (ABS0366)
PRODUCT REFERENCE	ARO-A17634
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	TOR1A

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-TOR1A, Anti-Torsin-1A, Anti-DQ2