

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Thyroid hormone receptor beta 1 Antibody (N402/13)
PRODUCT REFERENCE	ARO-A17952
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Thyroid hormone receptor beta 1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Thyroid hormone receptor beta 1, Anti-THRB, Anti-ERBA2