

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-GluK5/Grik5/KA2 kainate receptor Antibody (N279A/79)
PRODUCT REFERENCE	ARO-A18210
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	KA2 kainate receptor

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-KA2 kainate receptor, Anti-GRIK5, Anti-GluK5