

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Kv9.1 K ⁺ channel Antibody (N461/19)
PRODUCT REFERENCE	ARO-A18399
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	Kv9.1 K ⁺ channel

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Kv9.1 K ⁺ channel, Anti-KCNS1