

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PDCoV S1/Spike glycoprotein 1 (CTD) Polyclonal Antibody
PRODUCT REFERENCE	ARO-A19972
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	S1/Spike glycoprotein 1 (CTD)

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-S1/Spike glycoprotein 1 (CTD)