

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	BiTE Format Isotype Control for Bispecific Antibody
PRODUCT REFERENCE	ARO-A20089
CLONALITY	Monoclonal Antibody
ISOTYPE	VH-VL-VH'-VL'
TARGET	Bispecific Antibody

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Bispecific Antibody, Anti-Isotype Control, Anti-BiTE, T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain, Lysozyme C, 3.2.1.17, 1,4-beta-N-acetylmuramidase C, Allergen Gal d IV, Gal d 4, LYZ