

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Human IgG1 (LALA P329G) $\hat{\Delta}$ $\hat{\Delta}$ Isotype Control Antibody (13R4)
<b>PRODUCT REFERENCE</b>	ARO-A20528
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1 LALA P329G
<b>TARGET</b>	$\hat{\Gamma}^2$ -Galactosidase

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Mus musculus (Mouse)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti- $\hat{\Gamma}^2$ -Galactosidase, LALA PG,LALA,PG,LALA-PG,L234A, L235A, P329G (LALA-PG)