

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human GBA NANOBODY [®] , protein (SAA2277)
PRODUCT REFERENCE	ARO-A21108
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	GBA

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Lysosomal acid glucosylceramidase, Lysosomal acid GCCase, EC:3.2.1.45, Acid beta-glucosidase, Alglucerase, Beta-glucocerebrosidase, Beta-GC, Beta-glucosylceramidase 1, Cholesterol glucosyltransferase, SGTase, EC:2.4.1.-, Cholesteryl-beta-glucosidase,