

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Influenza A Virus H7 subtype Hemagglutinin/HA Antibody (P52E03)
PRODUCT REFERENCE	ARO-A21113
CLONALITY	Monoclonal Antibody
ISOTYPE	
TARGET	Influenza A Virus H7 subtype Hemagglutinin

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Hemagglutinin (H7 subtype)