



# Proteogenix

**Proteogenix S.A.S**

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

**Proteogenix INC.**

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

**PRODUCT NAME** Anti-CDH3/MSLN BiTE Antibody (AMG 305)

**PRODUCT REFERENCE** ARO-A21159

**CLONALITY** Monoclonal Antibody

**ISOTYPE**

**TARGET** CDH3

### PRODUCT INFORMATION

---

**ORIGIN SPECIES** *Mus musculus (Mouse)*

**FORM** Liquid

**BUFFER**

**STABILITY & STORAGE** Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

**ALIAS** Anti-Cadherin-3, Placental cadherin, P-cadherin, CDH3, CDHP, Mesothelin, CAK1 antigen, Pre-pro-megakaryocyte-potentiating factor, Megakaryocyte-potentiating factor, MPF, Mesothelin cleaved form, MSLN, MPF

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**