



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME Anti-FHV-1 glycoprotein C/gC Polyclonal Antibody

PRODUCT REFERENCE ARO-A21329

CLONALITY Polyclonal Antibody

ISOTYPE

TARGET envelope glycoprotein C, UL44

PRODUCT INFORMATION

ORIGIN SPECIES *Oryctolagus cuniculus (Rabbit)*

FORM Liquid

BUFFER

STABILITY & STORAGE Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS Anti-FHV-1 glycoprotein C

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY