

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-POLR2E Polyclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A21390
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	
<b>TARGET</b>	DNA-directed RNA polymerases I, II, and III subunit RPABC1, RNA polymerases I, II, and III subunit ABC1, DNA-directed RNA polymerase II 23 kDa polypeptide, DNA-directed RNA polymerase II subunit E, RPB5 homolog, XAP4, POLR2E

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-POLR2E