

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-ADH7 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21415
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	All-trans-retinol dehydrogenase [NAD(+)] ADH7, 1.1.1.105, Alcohol dehydrogenase class 4 mu/sigma chain, 1.1.1.1, Alcohol dehydrogenase class IV mu/sigma chain, Gastric alcohol dehydrogenase, Omega-hydroxydecanoate dehydrogenase ADH7, 1.1.1.66, Retinol deh

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ADH7