

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-POLR2K Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21480
CLONALITY	Polyclonal Antibody
ISOTYPE	
TARGET	DNA-directed RNA polymerases I, II, and III subunit RPABC4, RNA polymerases I, II, and III subunit ABC4, ABC10-alpha, DNA-directed RNA polymerase II subunit K, RNA polymerase II 7.0 kDa subunit, RPB7.0, RPB10alpha, POLR2K

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-POLR2K