



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Anti-NTMT2 Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A21606 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | |
| TARGET | N-terminal Xaa-Pro-Lys N-methyltransferase 2, 2.1.1.299, Alpha N-terminal protein methyltransferase 1B, Methyltransferase-like protein 11B, X-Pro-Lys N-terminal protein methyltransferase 1B, NTM1B, NTMT2, C1orf184, METTL11B, NRMT2 |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| ORIGIN SPECIES | <i>Oryctolagus cuniculus (Rabbit)</i> |
| FORM | Liquid |
| BUFFER | |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-NTMT2 |

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY