

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human MS4A2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21719
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	High affinity immunoglobulin epsilon receptor subunit beta, FcERI, Fc epsilon receptor I beta-chain, IgE Fc receptor subunit beta, Membrane-spanning 4-domains subfamily A member 2, MS4A2, APY, FCER1B, IGER

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-MS4A2