

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-FHIT Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21724
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Bis(5'-adenosyl)-triphosphatase, 3.6.1.29, AP3A hydrolase, AP3Aase, Adenosine 5'-monophosphoramidase FHIT, 3.9.1.-, Adenylylsulfatase, 3.6.2.1, Adenylylsulfate-ammonia adenylyltransferase, 2.7.7.51, Diadenosine 5',5'''-P1,P3-triphosphate hydrolase, Dinucl

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-FHIT