

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human NMRK1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21742
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Nicotinamide riboside kinase 1, NRK 1, NmR-K 1, 2.7.1.22, Nicotinic acid riboside kinase 1, 2.7.1.173, Ribosylnicotinamide kinase 1, RNK 1, Ribosylnicotinic acid kinase 1, NMRK1, C9orf95, NRK1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-NMRK1