

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-UBE2R2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21743
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Ubiquitin-conjugating enzyme E2 R2, 2.3.2.23, E2 ubiquitin-conjugating enzyme R2, Ubiquitin carrier protein R2, Ubiquitin-conjugating enzyme E2-CDC34B, Ubiquitin-protein ligase R2, UBE2R2, CDC34B, UBC3B

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-UBE2R2