

## PRODUCT DATASHEET

### OVERVIEW

---

|                          |                                                                                                                                                              |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>PRODUCT NAME</b>      | Anti-SMU1 Polyclonal Antibody                                                                                                                                |
| <b>PRODUCT REFERENCE</b> | ARO-A21745                                                                                                                                                   |
| <b>CLONALITY</b>         | Polyclonal Antibody                                                                                                                                          |
| <b>ISOTYPE</b>           | IgG                                                                                                                                                          |
| <b>TARGET</b>            | WD40 repeat-containing protein SMU1, Smu-1 suppressor of mec-8 and unc-52 protein homolog, WD40 repeat-containing protein SMU1, N-terminally processed, SMU1 |

### PRODUCT INFORMATION

---

|                                |                                                                                                                                                                                                                                                                                                            |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ORIGIN SPECIES</b>          | <i>Oryctolagus cuniculus (Rabbit)</i>                                                                                                                                                                                                                                                                      |
| <b>FORM</b>                    | Liquid                                                                                                                                                                                                                                                                                                     |
| <b>BUFFER</b>                  |                                                                                                                                                                                                                                                                                                            |
| <b>STABILITY &amp; STORAGE</b> | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| <b>ALIAS</b>                   | Anti-SMU1                                                                                                                                                                                                                                                                                                  |