

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-TSN Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21748
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Translin, 3.1.-.-, Component 3 of promoter of RISC, C3PO, TSN

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-TSN