

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-GTF2H2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21753
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	General transcription factor IIH subunit 2, Basic transcription factor 2 44 kDa subunit, BTF2 p44, General transcription factor IIH polypeptide 2, TFIIH basal transcription factor complex p44 subunit, GTF2H2, BTF2P44

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-GTF2H2