

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-PACC1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21791
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Proton-activated chloride channel, PAC, hPAC, Acid-sensitive outwardly-rectifying anion channel, ASOR, Proton-activated outwardly rectifying anion channel, PAORAC, Transmembrane protein 206, hTMEM206, PACC1, C1orf75, TMEM206

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-PACC1