



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-ST3GAL4 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21795
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 4, Alpha 2,3-ST 4, Beta-galactoside alpha-2,3-sialyltransferase 4, 2.4.3.2, 2.4.3.4, Alpha 2,3-sialyltransferase IV, Gal-NAc6S, Gal-beta-1,3-GalNAc-alpha-2,3-sialyltransferase, Gal-bet

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-ST3GAL4

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY