

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human IL1F10 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21811
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Interleukin-1 family member 10, IL-1F10, Family of interleukin 1-theta, FIL1 theta, Interleukin-1 HY2, IL-1HY2, Interleukin-1 theta, IL-1 theta, Interleukin-38, IL-38, IL1F10, FIL1T, IL1HY2, IL38, FKSG75, UNQ6119/PRO20041

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-IL1F10