

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-MEX3C Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21837
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	RNA-binding E3 ubiquitin-protein ligase MEX3C, 2.3.2.27, RING finger and KH domain-containing protein 2, RING finger protein 194, RING-type E3 ubiquitin transferase MEX3C, MEX3C, RKHD2, RNF194, BM-013

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-MEX3C