

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human SNAPC1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21842
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	snRNA-activating protein complex subunit 1, SNAPc subunit 1, Proximal sequence element-binding transcription factor subunit gamma, PSE-binding factor subunit gamma, PTF subunit gamma, Small nuclear RNA-activating complex polypeptide 1, snRNA-activating pr

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-SNAPC1