

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-MAN2A1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21843
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Alpha-mannosidase 2, 3.2.1.114, Golgi alpha-mannosidase II, AMan II, Man II, Mannosidase alpha class 2A member 1, Mannosyl-oligosaccharide 1,3-1,6-alpha-mannosidase, MAN2A1, MANA2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Oryctolagus cuniculus (Rabbit)</i>
FORM	Liquid
BUFFER	
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-MAN2A1