

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-Human ART4 Polyclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A21848
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	IgG
<b>TARGET</b>	Ecto-ADP-ribosyltransferase 4, 2.4.2.31, ADP-ribosyltransferase C2 and C3 toxin-like 4, ARTC4, Dombrock blood group carrier molecule, Mono(ADP-ribosyl)transferase 4, NAD(P)(+)-arginine ADP-ribosyltransferase 4, CD297, ART4, DO, DOK1

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Oryctolagus cuniculus (Rabbit)</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
<b>ALIAS</b>	Anti-ART4