



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Zebrafish HSP70 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21972
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	HSP70

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Danio rerio</i> , <i>Homo sapiens</i> , <i>Xenopus laevis</i> , <i>Brugia malayi</i> , <i>Cricetulus griseus</i> , <i>Pleurodeles waltl</i> , <i>Sus scrofa</i> , <i>Pongo abelii</i> , <i>Oryzias latipes</i> , <i>Onchocerca volvulus</i> , <i>Penicillium citrinum</i> , <i>Bos taurus</i> , <i>Gallus gallus</i> , <i>Oncorhynchus tshawytscha</i> , <i>Macaca fascicularis</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-stress protein HSP70, HSP70

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY