

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Zebrafish Amyloid-beta A4 protein/app Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21976
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Amyloid-beta, Amyloid-beta precursor protein, APP, ABPP, APPI, Alzheimer disease amyloid A4 protein homolog, Alzheimer disease amyloid protein, Amyloid precursor protein, Amyloid-beta (A4) precursor protein, Amyloid-beta A4 protein, Cerebral vascular amyl

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Danio rerio</i> , <i>Tetraodon fluviatilis</i> , <i>Takifugu rubripes</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-amyloid beta (A4) precursor protein a precursor, appa