

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human XIAP/BIRC3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21980
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	XIAP/BIRC3

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Rattus norvegicus, Mus musculus, Homo sapiens</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-E3 ubiquitin-protein ligase XIAP, 2.3.2.27, Baculoviral IAP repeat-containing protein 4, IAP-like protein, ILP, hILP, Inhibitor of apoptosis protein 3, IAP-3, hIAP-3, hIAP3, RING-type E3 ubiquitin transferase XIAP, X-linked inhibitor of apoptosis pro