

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human GPX2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A21997
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	GPX2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens, Mus musculus, Rattus norvegicus, Callithrix jacchus, Hylobates lar, Macaca fuscata fuscata, Pongo pygmaeus, Sapajus apella</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Glutathione peroxidase 2, GPx-2, GSHPx-2, EC:1.11.1.9, Gastrointestinal glutathione peroxidase, Glutathione peroxidase-gastrointestinal, GPx-GI, GSHPx-GI, Glutathione peroxidase-related protein 2, GPRP-2, Phospholipid hydroperoxide glutathione peroxi