

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human NSMCE2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22008
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	NSMCE2

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens, Bos taurus</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-E3 SUMO-protein ligase NSE2, EC:2.3.2.-, E3 SUMO-protein transferase NSE2, MMS21 homolog, hMMS21, Non-structural maintenance of chromosomes element 2 homolog, Non-SMC element 2 homolog, NSMCE2