

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human BAAT Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22014
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	BAAT

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Bile acid-CoA:amino acid N-acyltransferase, BACAT, BAT, EC:2.3.1.65, Bile acid-CoA thioesterase, Choloyl-CoA hydrolase, EC:3.1.2.27, Glycine N-choloyltransferase, Long-chain fatty-acyl-CoA hydrolase, EC:3.1.2.2, BAAT