

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-ALG8 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22016
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	ALG8

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens, Mus musculus, Bos taurus</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Dolichyl pyrophosphate Glc1Man9GlcNAc2 alpha-13-glucosyltransferase, EC:2.4.1.265, Asparagine-linked glycosylation protein 8 homolog, Dol-P-Glc:Glc(1)Man(9)GlcNAc(2)-PP-dolichyl alpha-13-glucosyltransferase, Dolichyl-P-Glc:Glc1Man9GlcNAc2-PP-dolichyl