

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CIART Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22017
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	CIART

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human, Mouse</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.
ALIAS	Anti-Circadian-associated transcriptional repressor, ChIP-derived repressor of network oscillator, Chrono, Computationally highlighted repressor of the network oscillator, CIART