

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|-------------------------------------|
| PRODUCT NAME | Anti-LCMV Glycoprotein GP2 Antibody |
| PRODUCT REFERENCE | ARO-A22032 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | GP2 |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Lymphocytic choriomeningitis virus</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |
| ALIAS | Anti-Pre-glycoprotein polyprotein GP complex, Pre-GP-C, Stable signal peptide, SSP, Glycoprotein G1, GP1, Glycoprotein G2, GP2, GPC |