



# ProteoGenix

## ProteoGenix S.A.S

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## ProteoGenix INC.

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Anti-CRHR1 Polyclonal Antibody
<b>PRODUCT REFERENCE</b>	ARO-A22051
<b>CLONALITY</b>	Polyclonal Antibody
<b>ISOTYPE</b>	IgG
<b>TARGET</b>	Corticotropin-releasing factor receptor 1; CRF-R-1; CRF-R1; CRFR-1; Corticotropin-releasing hormone receptor 1; CRH-R-1; CRH-R1; CRHR1; CRFR; CRFR1; CRHR

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	<i>Human, Mouse, Rat, Sheep</i>
<b>FORM</b>	Liquid
<b>BUFFER</b>	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
<b>STABILITY &amp; STORAGE</b>	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

### ALIAS

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**