

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Human CCL3L1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22067
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	C-C motif chemokine 3-like 1; G0/G1 switch regulatory protein 19-2; LD78-beta(1-70); PAT 464.2; Small-inducible cytokine A3-like 1; Tonsillar lymphocyte LD78 beta protein; LD78-beta(3-70); LD78-beta(5-70); CCL3L1; D17S1718; G0S19-2; SCYA3L1; CCL3L3

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Homo sapiens (Human)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS