



ProteoGenix

ProteoGenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

ProteoGenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Anti-Human NDUFB6 Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A22070 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6; Complex I-B17; CI-B17; NADH-ubiquinone oxidoreductase B17 subunit; NDUFB6 |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| ORIGIN SPECIES | <i>Homo sapiens (Human)</i> |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY