



Proteogenix

Proteogenix S.A.S

19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

products@proteogenix.fr

Proteogenix INC.

1000 Bearcat Way, Suite
100, Morrisville, NC 27560,
USA

Tel +1(919)234-1277

products@proteogenix.fr

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-CD197/CCR7 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22109
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	C-C chemokine receptor type 7; C-C CKR-7; CC-CKR-7; CCR-7; Epstein-Barr virus-induced G-protein coupled receptor 1; EB11; EBV-induced G-protein coupled receptor 1; MIP-3 beta receptor; CD197; Ccr7; Cmkbr7; Ebi1; Ebi1h

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mouse, Bovine, Human</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS

You can contact us at products@proteogenix.fr for any question

OUR PRODUCTS ARE FOR RESEARCH USE ONLY