

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse CD121a/IL1R1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22118
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Interleukin-1 receptor type 1; IL-1R-1; IL-1RT-1; IL-1RT1; 3.2.2.6; CD121 antigen-like family member A; Interleukin-1 receptor alpha; IL-1R-alpha; Interleukin-1 receptor type I; p80; CD121a; Interleukin-1 receptor type 1, membrane form; mL-1R1; mL-1RI;

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS