

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-Mouse RIPK3 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22153
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Receptor-interacting serine/threonine-protein kinase 3; 2.7.11.1; RIP-like protein kinase 3; Receptor-interacting protein 3; RIP-3; mRIP3; Ripk3; Rip3

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Mus musculus (Mouse)</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS