

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-TMSB4X Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22166
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Thymosin beta-4; T beta 4; Hemoregulatory peptide AcSDKP; N-acetyl-SDKP; AcSDKP; Seraspenide; Tmsb4x; Ptmb4; Tmsb4

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Pig, Rat, Frog, Human, Horse, Mouse, Rabbit, Bovine</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS