

PRODUCT DATASHEET

OVERVIEW

| | |
|--------------------------|--|
| PRODUCT NAME | Anti-Der f 2 Polyclonal Antibody |
| PRODUCT REFERENCE | ARO-A22211 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | Mite group 2 allergen Der f 2, Allergen Der f II, Der f 2, DERF2 |

PRODUCT INFORMATION

| | |
|--------------------------------|---|
| ORIGIN SPECIES | <i>Dermatophagoides farinae</i> (American house dust mite) |
| FORM | Liquid |
| BUFFER | 0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300 |
| STABILITY & STORAGE | Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs. |

ALIAS