

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-FPV NS1 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22274
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Initiator protein NS1, NS1, EC:3.1.21.-, EC:3.6.4.12, NCVP1, Non-capsid protein NS-1, Non-structural protein 1, Non-structural protein NS1, NS1

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Feline panleukopenia virus</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS