

PRODUCT DATASHEET

OVERVIEW

PRODUCT NAME	Anti-VP2 Polyclonal Antibody
PRODUCT REFERENCE	ARO-A22312
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	1B VP2 mature peptide; 1B VP2 mature peptide

PRODUCT INFORMATION

ORIGIN SPECIES	<i>Human hepatitis A virus genotype IIB, Human hepatitis A virus genotype IIIB, Human hepatitis A virus genotype IIA, Human hepatitis A virus genotype IIIA, Human hepatitis A virus genotype IA, Simian hepatitis A virus genotype IV, Simian hepatitis A virus g</i>
FORM	Liquid
BUFFER	0.01M PBS (pH 7.4), 50% Glycerol, 0.05% Proclin 300
STABILITY & STORAGE	Store at 2°C to 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

ALIAS